**PoloDent Instrument Disinfection** is a first rate aldehyde- and phenol-free concentrate for the cleaning and disinfection of general dental instruments (mirrors, forceps etc.), as well as alkali- and alcohol-sensitive rotary instruments and endodontic instruments. This product is 99% biodegradable.

**Composition:** Contains amongst other ingredients: Less than 5% Amphoteric surfactant, Disinfectants, Parfum, Chlorhexidine Digluconate, Didecylmethyl ammonium chloride, Benzalkonium chloride.

**Effective spectrum:** Bactericidal, fungicidal, tuberculocidal, virucidal, E-coli, HIV1, Human Influenza B Virus (HIBV), MRSA, gram negative- and gram positive bacteria. Clostridium Difficile vegetative cell (growing cells) of Gram+ organisms. Tested to: EN 13624, EN 13727, EN 14476, EN 14204, EN 13704.

- Bactericidal in 60 minutes (2% - 20 ml)
- Fungicidal in 60 minutes (2% - 20 ml)
- Tuberculocidal in 60 minutes (2% - 20 ml) and 30 minutes (4% - 40 ml)
- Virucidal in 15 minutes (5% - 50 ml)
- HIV1/HIBV in 15 minutes (2% - 20 ml)
- MRSA in 30 seconds (2% - 20 ml)

**Directions for use:** Use as a 2% or 4% solution. For use in ultrasonic equipment, we recommend a 4% solution. Do not mix with ordinary household detergents. Make up with cold or luke warm water. Immerse objects completely and ensure that all surfaces are wet and that the air is expelled. After the required contact time rinse objects thoroughly under running water (or sterile water). Dry and/or sterilize instruments immediately. The maximum standing time for the ready-for-use solution is 7 days. In case of high contamination with blood, saliva etc. we recommend renewal of the solution 2-3 times per week.

**Caution:** R34: Causes burns. S2: Keep out of reach of children. S25: Avoid contact with eyes. S26: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. S28: After contact with skin, wash immediately with plenty of water. S29: Do not empty into drains.

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Polo MB Oisterwijk B.V.